

HCG antibody (orb538856)

Western blot analysis of Cx36 expression in various tissues. The blot shows protein bands for Cx36 (indicated by an arrowhead) in U-87MG, U-251MG, MCF7, Mouse testis, Rat ovary, and Rat testis. Molecular weight markers are indicated on the left (70kDa, 55kDa, 40kDa, 35kDa, 25kDa, 15kDa, 10kDa). The Cx36 band is prominent in the U-87MG, U-251MG, and MCF7 lanes, and also visible in the Mouse testis, Rat ovary, and Rat testis lanes.

Western blot
analysis of
extracts of
var...