


Product Datasheet

CHRNA3 antibody (orb538893)

Description	Rabbit polyclonal antibody to CHRNA3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CHRNA3 (NP_000740.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	55kDa
Uniprot ID	Q05901
Entrez	1142
Dilution Range	WB: 1:1000-1:4000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of various tissues.