

Product Datasheet

Acss2 antibody (orb53891)

Description	Acss2 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Mouse Acetyl-coenzyme A synthetase, cytoplasmic protein (1-701AA)
Target	Acss2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9QXG4
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.