



Product Datasheet

CLCA1 antibody (orb538916)


Description	Rabbit polyclonal antibody to CLCA1
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Recombinant fusion protein of human CLCA1 (NP_001276.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	85kDa
Uniprot ID	A8K7I4
Entrez	1179
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse col...



Immunofluorescence analysis of human col...



Immunofluorescence analysis of rat colon...