



Product Datasheet

CMAS antibody (orb538948)


Description	Rabbit polyclonal antibody to CMAS
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CMAS(NP_061156.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48kDa
Uniprot ID	Q8NFW8
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse brain...



Immunohistochemical staining of rat brain...



Western blot analysis of extracts of various...