



Product Datasheet

COPS5 antibody (orb538980)


Description	Rabbit polyclonal antibody to COPS5
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human JAB1/CSN5/JAB1/CSN5/COPS5(NP_006828.2).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37kDa
Uniprot ID	Q92905
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of NIH/3T3 cells



Immunofluorescence analysis of C6 cells



Immunofluorescence analysis of HeLa cells