



Product Datasheet

CSF3 antibody (orb539037)


Description	Rabbit polyclonal antibody to CSF3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human CSF3(NP_000750.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22kDa
Uniprot ID	P09919
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse spleen



Immunohistochemical staining of human colon



Immunohistochemical staining of rat spleen