



Product Datasheet

CTGF antibody (orb539055)


Description	Rabbit polyclonal antibody to CTGF
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human CTGF(NP_001892.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35KDa
Uniprot ID	P29279
Dilution Range	WB: 1:200-1:500, IHC-P: 1:50-1:200, IF: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of C6 cells ...



Immunohistochemical staining of mouse ki...



Immunohistochemical staining of human br...