



Product Datasheet

CTLA4 antibody (orb539056)

Description	Rabbit polyclonal antibody to CTLA4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CD152/CTLA-4(NP_005205.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36kDa
Uniprot ID	P16410
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:10-1:100, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human to...



Western blot analysis of extracts of var...