



Product Datasheet

CTSA antibody (orb539063)

Description	Rabbit polyclonal antibody to CTSA
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CTSA(NP_001121167.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54KDa
Uniprot ID	P10619
Dilution Range	WB: 1:500-1:2000, IF: 1:20-1:50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cells.



Western blot analysis of extracts of various tissues.