


Product Datasheet

CWC22 antibody (orb539077)

Description	Rabbit polyclonal antibody to CWC22
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CWC22(NP_065994.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110kDa
Uniprot ID	Q9HCG8
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...