



Product Datasheet

CYP11B2 antibody (orb539089)


Description	Rabbit polyclonal antibody to CYP11B2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human CYP11B2(NP_000489.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Uniprot ID	P19099
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse ad...



Immunohistochemical staining of rat adre...



Western blot analysis of extracts of var...