



Product Datasheet

CYP20A1 antibody (orb539091)

Description	Rabbit polyclonal antibody to CYP20A1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CYP20A1 (NP_803882.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	46kDa
Uniprot ID	Q6UW02
Entrez	57404
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat kidn...



Western blot analysis of extracts of var...