

Product Datasheet

CYP2A6 antibody (orb539095)

Description	Rabbit polyclonal antibody to CYP2A6
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CYP2A6(NP_000753.3).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50kDa
Uniprot ID	P11509
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.