


Product Datasheet

DDX46 antibody (orb539148)

Description	Rabbit polyclonal antibody to DDX46
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human DDX46(NP_055644.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	117kDa
Uniprot ID	Q7L014
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...