


## Product Datasheet

### DDX58 antibody (orb539152)

|                            |                                                               |
|----------------------------|---------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to DDX58                           |
| <b>Reactivity</b>          | Human, Mouse, Rat                                             |
| <b>Conjugation</b>         | Unconjugated                                                  |
| <b>Tested Applications</b> | ICC, IF, WB                                                   |
| <b>Immunogen</b>           | Recombinant fusion protein of human RIG-I/DDX58(NP_055129.2). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.            |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.                   |
| <b>Note</b>                | For research use only                                         |
| <b>Isotype</b>             | IgG                                                           |
| <b>Clonality</b>           | Polyclonal                                                    |
| <b>MW</b>                  | 102-110kDa                                                    |
| <b>Uniprot ID</b>          | <b>095786</b>                                                 |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                               |



Western blot analysis of extracts of var...