



Product Datasheet

DDX6 antibody (orb539153)


Description	Rabbit polyclonal antibody to DDX6
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human DDX6(NP_001244120.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54kDa
Uniprot ID	P26196
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human mammary gland tissue.



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of human lung tissue.