


Product Datasheet

DNAJC24 antibody (orb539215)

Description	Rabbit polyclonal antibody to DNAJC24
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human DNAJC24(NP_859057.4).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	17kDa
Uniprot ID	Q6P3W2
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...