


## Product Datasheet

### DNAL1 antibody (orb539217)

|                            |                                                         |
|----------------------------|---------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to DNAL1                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                       |
| <b>Conjugation</b>         | Unconjugated                                            |
| <b>Tested Applications</b> | IHC, WB                                                 |
| <b>Immunogen</b>           | Recombinant fusion protein of human DNAL1(NP_113615.2). |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.      |
| <b>Storage</b>             | Store at -20°C. Avoid freeze / thaw cycles.             |
| <b>Note</b>                | For research use only                                   |
| <b>Isotype</b>             | IgG                                                     |
| <b>Clonality</b>           | Polyclonal                                              |
| <b>MW</b>                  | 22kDa                                                   |
| <b>Uniprot ID</b>          | <a href="#">Q4LDG9</a>                                  |
| <b>Dilution Range</b>      | WB: 1:500-1:2000, IHC-P: 1:50-1:200                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                         |



Western blot analysis of extracts of var...