



Product Datasheet

DNM1L antibody (orb539220)

Description	Rabbit polyclonal antibody to DNM1L
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human DNM1L(NP_036192.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82kDa
Uniprot ID	000429
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of C6 cells

...



Western blot analysis of extracts of var...