



Product Datasheet

DNMT3A antibody (orb539226)

Description	Rabbit polyclonal antibody to DNMT3A
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human DNMT3A(NP_072046.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	140kDa
Uniprot ID	Q9Y6K1
Dilution Range	WB: 1:500-1:2000, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Western blot analysis of extracts of var...