


Product Datasheet

B4GALT7 antibody (orb53923)

Description	B4GALT7 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Beta-1,4-galactosyltransferase 7 protein (52-327AA)
Target	B4GALT7
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBV7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse brain tis...