


Product Datasheet

DPYSL5 antibody (orb539247)

Description	Rabbit polyclonal antibody to DPYSL5
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human CRMP5/CRMP5/DPYSL5 (NP_064519.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	61kDa
Uniprot ID	Q9BPU6
Entrez	56896
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...