

Product Datasheet

DSTYK antibody (orb539254)

Description	Rabbit polyclonal antibody to DSTYK
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human DSTYK(NP_955749.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105kDa
Uniprot ID	Q6XUX3
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.