



Product Datasheet

DUT antibody (orb539266)


Description	Rabbit polyclonal antibody to DUT
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human DUT(NP_001020419.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24kDa
Uniprot ID	P33316
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells



Immunohistochemical staining of human appendix



Immunohistochemical staining of human liver