



Product Datasheet

EGLN1/EGLN2 antibody (orb539302)


Description	Rabbit polyclonal antibody to EGLN1/EGLN2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human EGLN1/EGLN2(NP_071334.1/NP_444274.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	44kDa
Uniprot ID	Q9GZT9, Q96KS0
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human gastric tissue.



Immunohistochemical staining of human stomach tissue.



Immunohistochemical staining of rat heart tissue.