


Product Datasheet

EGR1 antibody (orb539304)

Description	Rabbit polyclonal antibody to EGR1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human EGR1 (NP_001955.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	75-80KDa
Uniprot ID	P18146
Entrez	1958
Dilution Range	WB: 1:200-1:1000, IHC-P: 1:20-1:200, IF: 1:10-1:100, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...