



Product Datasheet

EIF3H antibody (orb539315)

Description	Rabbit polyclonal antibody to EIF3H
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Recombinant fusion protein of human EIF3H(NP_003747.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa
Uniprot ID	O15372
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF7 cell...



Western blot analysis of extracts of var...