



Product Datasheet

EZH2 antibody (orb539409)


Description	Rabbit polyclonal antibody to EZH2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human EZH2/KMT6 (NP_001190176.1).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	98KDa
Uniprot ID	Q15910
Entrez	2146
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse li...



Immunohistochemical staining of rat lung...



Western blot analysis of extracts of var...