



Product Datasheet

Tissue factor antibody (orb539415)


Description	Rabbit Polyclonal antibody to Tissue factor
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Coagulation factor III/Tissue Factor(NP_001984.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53KDa/45KDa
Uniprot ID	P13726
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human esophagus tissue.



Immunohistochemical staining of rat liver tissue.



Western blot analysis of extracts of various cell lines.