


Product Datasheet

FAM107A antibody (orb539428)

Description	Rabbit polyclonal antibody to FAM107A
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human FAM107A(NP_009108.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	18kDa
Uniprot ID	O95990
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...