


Product Datasheet


FBL antibody (orb539447)

Description	Rabbit polyclonal antibody to FBL
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Recombinant fusion protein of human Fibrillarin/U3 RNP (NP_001427.2).
Preservatives	Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purification
MW	37kDa
Uniprot ID	P22087
Entrez	2091
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, IP: 1:20-1:50
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of C6 cells

...



Immunofluorescence analysis of U2OS cell...



Immunofluorescence analysis of HeLa cell...