



Product Datasheet

FGF12 antibody (orb539478)


Description	Rabbit polyclonal antibody to FGF12
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human FGF12(NP_004104.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20kDa
Uniprot ID	P61328
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse testis tissue.



Immunohistochemical staining of human colon tissue.



Immunohistochemical staining of rat testis tissue.