



Product Datasheet

FLCN antibody (orb539494)

Description	Rabbit polyclonal antibody to FLCN
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human FLCN(NP_659434.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70KDa
Uniprot ID	Q8NFG4
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of A549 cell...



Western blot analysis of extracts of var...