

Product Datasheet

FOXP1 antibody (orb539524)

Description	Rabbit polyclonal antibody to FOXP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human FOXP1(NP_116071.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	75KDa
Uniprot ID	Q9H334
Dilution Range	WB: 1:1000-1:3000
Expiration Date	12 months from date of receipt.

