



Product Datasheet

FUT4 antibody (orb539544)


Description	Rabbit polyclonal antibody to FUT4
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human FUT4(NP_002024.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59KDa
Uniprot ID	P22083
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat pancreas.



Immunohistochemical staining of human liver.



Immunohistochemical staining of human colon.