


Product Datasheet

GABBR1 antibody (orb539567)

Description	Rabbit polyclonal antibody to GABBR1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human GABBR1(NP_068703.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	95kDa
Uniprot ID	Q9UBS5
Dilution Range	WB: 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...