



Product Datasheet

GCC2 antibody (orb539597)


Description	Rabbit polyclonal antibody to GCC2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human GCC2(NP_852118.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	196kDa
Uniprot ID	Q8IWJ2
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human ap...



Immunohistochemical staining of human st...



Immunohistochemical staining of human co...