



Product Datasheet

GLP1 antibody (orb539598)


Description	Rabbit Polyclonal antibody to GLP1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	Recombinant fusion protein of human Pro-glucagon(NP_002045.1).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21kDa
Uniprot ID	P01275
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse pancreas tissue



Immunofluorescence analysis of rat pancreas tissue



Immunohistochemical staining of mouse pancreas tissue