



Product Datasheet

GNAI3 antibody (orb539654)

Description	Rabbit polyclonal antibody to GNAI3
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, IP, WB
Immunogen	Recombinant fusion protein of human GNAI3(NP_006487.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45kDa
Uniprot ID	P08754
Dilution Range	WB: 1:500-1:2000, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from n...



Western blot analysis of extracts of var...