



Product Datasheet

GRIN1 antibody (orb539722)


Description	Rabbit polyclonal antibody to GRIN1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide of human NMDAR1(NP_015566.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105kDa
Uniprot ID	Q05586
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of SKOV3 cells...



Immunofluorescence analysis of rat brain...



Western blot analysis of extracts of various cell lines...