



Product Datasheet

GRM4 antibody (orb539729)


Description	Rabbit polyclonal antibody to GRM4
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human GRM4(NP_000832.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	115kDa
Uniprot ID	Q14833
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF-7 cells.



Immunohistochemical staining of human brain tissue.



Immunohistochemical staining of rat brain tissue.