


Product Datasheet

GSK3B antibody (orb539735)

Description	Rabbit polyclonal antibody to GSK3B
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human GSK3 β (NP_001139628.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48KDa
Uniprot ID	P49841
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:100, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...