



Product Datasheet

GSS antibody (orb539738)

Description	Rabbit polyclonal antibody to GSS
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human Glutathione Synthetase (Glutathione Synthetase (GSS))(NP_000169.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52KDa
Uniprot ID	P48637
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from n...



Western blot analysis of extracts of var...