


Product Datasheet

GYPC antibody (orb539760)

Description	Rabbit polyclonal antibody to GYPC
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Glycophorin C (Glycophorin C (GYPC))(NP_002092.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34kDa
Uniprot ID	P04921
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...