



Product Datasheet

HACE1 antibody (orb539769)


Description	Rabbit polyclonal antibody to HACE1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthetic peptide of human HACE1(NP_065822.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	102kDa
Uniprot ID	Q8IYU2
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat testis...



Immunohistochemical staining of human pl...



Immunohistochemical staining of human li...