



Product Datasheet

HAVCR2 antibody (orb539783)

Description	Rabbit polyclonal antibody to HAVCR2
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	A synthetic peptide of human TIM-3/HAVCR2(NP_116171.3).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45-70kDa
Uniprot ID	Q8TDQ0
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of Hela cell...



Western blot analysis of extracts of var...