



Product Datasheet

HBA1 antibody (orb539784)

Description	Rabbit polyclonal antibody to HBA1
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant fusion protein of human Hemoglobin subunit alpha (Hemoglobin subunit alpha (HBA1))(NP_000508.1).
Preservatives	PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa/25kDa
Uniprot ID	P69905
Dilution Range	WB: 1:500-1:2000, IF: 1:50-1:100
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF7 cell...



Western blot analysis of extracts of var...