


Product Datasheet

HDAC3 antibody (orb539796)

Description	Rabbit polyclonal antibody to HDAC3
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human HDAC3(O15379).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49kDa
Uniprot ID	O15379
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200, IP: 1:50-1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of var...