



Product Datasheet

HIGD1A antibody (orb539815)

Description	Rabbit polyclonal antibody to HIGD1A
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide of human HIGD1A(NP_054775.2).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	14kDa
Uniprot ID	Q9Y241
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U-2 OS cells



Western blot analysis of extracts of various cell lines